



VEE OneLab compared to VEE Pro

- [Alliances](#)
- [Education](#)
- [Events](#)
- [Library](#)
- [News Releases](#)
- [Order Information](#)
- [Technical Support](#)

The E2123G VEE OneLab 6.0 for Win95/98/NT/2000 is designed for the individual engineer or scientist doing PC measurements, data acquisition, or design.

The E2120G VEE Pro 6.0 for Win95/98/NT/2000 is designed to support more instruments at once, larger programs (and so includes Program Explorer, Profiler, and Find), and sharing programs among engineers (importing libraries of UserFunctions, and Merge).

VEE OneLab 6.0 provides a subset of the capability of VEE Pro 6.0. VEE OneLab includes all of the new VEE Pro features except the Callable VEE ActiveX Automation Server.

VEE OneLab 6.0 does not support the following VEE Pro 6.0 objects and features. VEE OneLab 6.0 cannot load a program containing any of these objects.

Unsupported Major Features

- Callable VEE ActiveX Automation Server
- Limit of 4 instruments configured in the Instrument Manager
- VEE Pro Runtime not included in the VEE OneLab license or bits
- Web Server

Unsupported Objects and Functions

- Sequencer object
- Import Library object (Import disabled for UserFunction and Remote Function; Import Compiled Function works)
- Execute Program object
- To/From DDE object
- To/From Socket object
- HP-UX objects (To/From Named Pipe, RMB/UX, stdin, etc.)

Unsupported Menu Picks and Features

- Merge menu pick
- Merge Library menu pick
- Save Objects menu pick
- Save Secured menu picks
- Object menu picks for Find, Locate, Calls, and Called By Find
- Program Explorer
- Profiler
- Call Stack

VEE OneLab is not available on the HP-UX platform.

